

## (USA) ENA HRi esthetic restorative system - Instructions for use

ENA HRi is a light curing radiopaque composite for direct and indirect esthetic restorations in anterior and posterior areas. It follows the standard ISO 4049:2000. The system includes:

### HRi Universal Enamel

To obtain an enamel composite material that acts like natural enamel it must be highly translucent and must have the same refractive index. When this index is matched, then thicker layers of composite enamel appears whiter (high value/high luminosity/low translucency) and when it is applied in a thinner layer it appears more translucent (low value/low luminosity/high translucency), while increasing the thickness of the enamel layers of standard composite increases the percentage of grey in proportion to its thickness (glass-like effect). The new **Universal Enamels (UE) have the same refractive index as the natural tooth** and the same high luminosity as natural enamel. These two unique properties require the material to be applied with a technique different from any other enamel shade of composite you may be familiar with. UE shades should be **applied at similar thickness, slightly thinner, as enamel on the tooth that is being restored**, showing no visible margin. **Thick layer of UE appear whiter (higher value), thin layer will be more transparent.** In the incisal area, where there is no dentin present, the enamel produces the blue-amber opalescent effect, as this multichromatic opalescent characteristics exhibited by natural enamel has been engineered into these Universal Enamel shades. **Note: If you want to further enhance the opalescence effect in the incisal area, you may use the Opalescent Blue Natural (OBN) or Opalescent Amber (OA)**, available separately; to reproduce white characterization area, use intensive white IM, IWS or IW, covering these bodies with a 0,3-0,5 mm layer of Universal Enamel UE (even thinner for enhancing the intensive) as thicker layers can cover these bodies. By carefully observing natural teeth it is possible to distinguish in the enamel different levels of translucency depending on the age of the patient. **Note: no matter which universal enamel is used, value can be increased by increasing thickness.** Three universal enamel shades are available:

<b>UE1</b>	low value in thin layer; with amber effects; increasing thickness the value increases
<b>UE2</b>	medium value that becomes high value increasing thickness
<b>UE3</b>	very high value, really white, to be used only for very white or bleached teeth

### Composition of Universal Enamel

- Monomer matrix: Diurethandimethacrylate, Iso-propyliden-bis (2(3)-hydroxy-3(2)-4(phenoxy)propyl)-bis(methacrylate)(Bis-GMA); 1,4 - Butandioldimethacrylate.
- Content of filler: 80% weight. Glass filler: mean particle size 1,0 µm, 12,5% weight. Nano zirconium oxide particles: particle size 20nm.

## Intensive

These shades are used for further characterization of the enamel (ridges and cusps) and are used within the surface of a Universal Enamel to mimic hypo calcification or other extremely white areas.

The intensive white shades represent demineralized enamel areas and can appear in all areas of the tooth (cervical, middle and incisal third).

<b>IM</b>	Intensive Milky	A Warm Opaque White
<b>IWS</b>	Intensive White Spot	An Intensive White intermediate
<b>IW</b>	Intensive White	A Cold Translucent White

## Dentin

A modern composite system has to include dentins with a fluorescence degree calibrated to the natural tooth.

The average **chroma** of natural teeth (central incisors, lateral incisors and canines) is in the region of 580 nm. The "A" shades of the Vita®\* shade guide are closer to the average chroma of natural teeth. For this reason, we developed the new Universal Dentin (UD) shades close to the **hue-chroma (chromaticity)** of natural teeth. These new Universal Dentin shades have a high brightness (higher value) and are calibrated to match the fluorescence and opacity of natural dentin. When determining the basic chroma of the tooth the areas that are most suitable are the cervical and middle third. In complex restorations, the final shade is created by using the basic hue and then next two darker dentins (for this UD5 and UD6 are available). In most restorations only one shade of dentin will be required, as margin is not visible thanks to the new universal enamels. New UD0 and UD0,5 shades are useful to restore very light or bleached teeth.

9 Fluorescent Dentin	UD0 - UD0,5 - UD1 (A1*) - UD2 (A2*) - UD3 (A3*) - UD3,5 (A3,5*) - UD4 (A4*) - UD5 - UD6
----------------------	---

## Composition of dentin and intensive

- Monomer matrix: Diurethandimethacrylate, Iso-propyliden-bis (2(3)-hydroxy-3(2)-4 (phenoxy) propyl)-bis(methacrylate)(Bis-GMA); 1,4 - Butandioldimethacrylate.

- Total content of fillers: 75% weight (53% volume); glass filler: mean particle size 0.7 µm; highly dispersed silicone dioxide: mean particle size 0.04 µm.

### Clinical indications

Class I (all cavities)	Class II (small and medium cavities)	Class III (all cavities)
Class IV (all cavities)	Class V (all cavities)	Sealing
Total and partial vestibular covering	Cosmetic corrections	Complex restorations
Inlays Class I (all cavities)	Inlays Class II (all cavities)	Inlays Class IV (all cavities)
Laminated veneers	Onlays	Restoration of prosthetic cores

### Contra-indications

Uncured resin could cause skin allergy: User should wear gloves. In case of known allergy to some of the components do not use it.

### Side effects

In deep cavities we suggest the use of a liner in order to avoid pulpal irritation.

### Materials to be avoided

Materials containing phenolics (like eugenol) could inhibit composite curing. Avoid the use of these materials as liners.

\* colors of Vita® shade guide. Vita® is a registered trademark of Vita Zahnfabrik H. Rauter mbH & Co. KG, Bad Säckingen - D

# DIRECT TECHNIQUE

FILLINGS AND DIRECT ESTHETIC RESTORATIONS OF CLASS I-II-III-IV-V.

## Preparation

- Clean with fluoride-free prophylaxis paste.
- Choose colors with Vita® shade guide or with ENA HRi composite shade guide, and fill in the "Color Chart".
- Preparation: for anterior teeth, use a conservative preparation with bevel, which allows a good enamel etching (for posterior restorations, no bevel is needed). We suggest ENA Shiny preparation kit of Dr. L. Vanini, where Shiny 33 rubber, for polishing preparation, is included.
- We suggest using a rubber dam.
- In case of interproximal restorations, use a transparent matrix.

## Etching and bonding

Follow your normal technique. We suggest 35%-38% phosphoric acid (ENA ETCH) for 15-30 seconds for enamel and vital dentin, 1 min. for sclerotic dentin and 2 min. for non-vital dentin. Wash and dry the etched surface with oil-free air; etched enamel looks white calcareous. Etched surfaces should not be contaminated before the application of bonding material (we recommend ENA BOND and Rock Bond, but ENA HRi works perfectly with your bonding system of choice). In case of contamination with saliva, wash, dry and etch again (avoid dehydrating the dentin). Apply a thin coat of bonding material on etched surfaces of dentin and enamel, pulling it down carefully on the margins and cure it for 40 sec. (with Translux CL or Nou-Lite halogen light curing units). (For ENA BOND apply and cure a second coat). Be careful not to contaminate the oxygen inhibition layer left after curing to assure a strong chemical bond to the composite.

## Composite application

Take ENA HRi out of the syringe or "tips"; we suggest to warm the composite in the ENA HEAT composite warmer to 39° C. Apply very small quantities of material by pulling it down with a brush (Micerium "M" brush for anterior and "F" for posteriors, and Micerium Silicone Brushes) in order to avoid any bubbles.

**NOTE: do not wet Universal Enamel with any resin or bonding because it will change the refraction index and cause the composite to become too opaque.**

Use a "waves" application technique in order to allow a better light diffusion effect. Follow the stratification technique described in the next paragraph. Cure layers of 1-1,5 mm (no more than 2 mm) for 40 seconds, from all sides of the build up; keep the light-curing tip as close as possible to the restoration. Oxygen leaves a thin layer of uncured composite: this layer should not be contaminated or wetted because it creates a chemical connection between the different layers of composite. We advise to apply an Air Block (Shiny G), when the restoration is finished and before the final light curing takes place. This glycerine based product eliminates the oxygen inhibition layer.

**CURING:** Working time under standard light is approximately 3 minutes. During a long procedure, cover the composite with an opaque foil or use a color palette with orange or black cover (COSSTAIN01). **NOTE:** avoid direct light of the overhead light and turn off if possible. Cure each layer for 40 seconds.


## Finishing and polishing

Use diamond burs and diamond pastes. Do not use any disc buccally in order to avoid destroying the texture surface. We suggest to use the complete finishing and polishing system ENA Shiny.

## Anatomic stratification technique of Dr. L. Vanini

In order to maximise the characteristics of the **ENA HRi** System, we suggest following the anatomic stratification technique of Dr. Lorenzo Vanini; any other stratification technique that does not respect the anatomy of a natural tooth would considerably limit the esthetic performance of this system.

### Anterior complex restorations



The lingual enamel is applied by means of a silicone matrix, using a Universal Enamel (UE1, UE2, UE3). The interproximal enamel wall is built up using the same Universal Enamel.






The thickness of all enamel layers is the same as if the natural enamel was still there.

To obtain a natural chromatic composition in complex restorations, two different dentin shades are used. Once the final shade is established, increase the dentin shade by two for the first layer. For example, if the desired shade is A1, the first shade used cervically would be UD3. This would be covered by UD1 to be applied more incisally, creating the mamelon shape.

Now, if desired, intensive shades (IW, IWS, IM) are built in and characterizations are reproduced with **ENA STAINS** (available in white, yellow, orange, blue, brown, dark brown) and opalescent amber (OA).

HRi gives a natural opalescent effect. If we need to reinforce the blue-amber opalescent effect, Opalescent shades OBN (Blue) and OA (Amber) can be used.

Finally the vestibular enamel is applied using a Universal Enamel.




### Complex Restorations (2 dentins, 1 enamel)




## Anterior medium and small restoration



Medium Size restorations  
(2 dentins, 1 enamel)




Common restorations  
(1 dentin, 1 enamel)




Enamel Only restorations  
(1 enamel)



## Posterior direct restoration



(1 or 2 dentin, 1 enamel)



## INDIRECT TECHNIQUE

INLAY, ONLAY AND VENEERS, COMBINATION PROSTHESIS AND IMPLANTS

**ENA HRi** can be used indirectly for anterior or posterior inlays, onlays, veneers, over implants and in combination cases. The dental technicians use **ENA HRi** with the same stratification technique as they use with modern ceramic systems.

### Preparation

Preparation should be made without undercuts, and for posterior restorations, slightly tapered diamonds are recommended to round out internal edges. Minimum thickness of composite layers should be >1.5 mm to avoid breakage.

Close undercuts by using **ENA HRi Flow** composite.

### Impression and Temporary

Take an impression and use **ENA Temp** for temporary inlay and cement it with eugenol-free cement. For inlays it is possible to use **ENA Soft** elastic composite. Its elastic properties allow for the complete and easy removal of the temporary inlay which leaves the preparation clean.

### Laboratory procedure

Pour a model with extra-hard plaster. After the plaster sets, remove the impression and apply an oil-free separator (**TEMP SEP**) to the model. Follow the same stratification technique as in the direct method. For inlays, first build up the external walls and then the occlusal areas. It is possible to use **ENA Stains** between Dentin and Universal Enamel. Each layer should not be thicker than 2 mm and should be cured for 40 sec.

Recommended final curing time is 11 minutes using a high power light curing unit such as **LaborluxL** or if using an 86W light box like **Lampadaplust** final curing time is 30 minutes. Finish with burs and polish with **ENA Shiny** brushes and diamond pastes. Wash with soap and water and dry with oil-free air spray.

**Note.** For further technical instructions also on restorations on metal and fiber structure, please consult the manual "ENA HRi Tender, laboratory procedures".

## Luting

Remove the temporary appliance and clean the preparation. Try-in the restoration carefully and proceed with any adjustments. Post-cure in an oven like Lampadaplust for 9 min. Apply the rubber dam. Clean the surface of the preparation with alcohol and sandblast. Etch the cavity and apply two coats ENA BOND but do not cure. Sandblast the internal part of the composite restoration, then clean it with alcohol; apply the bond resin but do not cure. Warm a small amount of ENA HRi enamel or a light dentin shade to 50°C in the warmer (ENA HEAT) and apply it to the inside of the restoration. When restoration is in place, apply a small amount of pressure either mechanically or manually. Remove composite excess at margins and cure for at least 80 seconds from each side of the tooth. Check the occlusion, finish and polish with ENA Shiny system, using burs, strips and diamond pastes.

**Note.** in case of inlay thickness over 2 mm use a dual luting composite such as ENA CEM (see instructions for details).

## Curing information

It is necessary to use a light-curing unit with a spectrum of 350 - 500 nm. The required physical results can be reached only if using multi-wall reflecting unit. For this reason we suggest a periodical check of the light intensity following the manufacturer's instructions. Most curing units reach a complete cure to a depth of 4.6 mm. Optimal values are reached at 2.3 mm.

Laboratory curing times:

- |  |                                    |
|--|------------------------------------|
| - LaborluxL (MICERIUM)                       | ca. 90 sec. (final curing 16 min.) |
| - Spektra 2000 (Schütz-Dental)               | ca. 90 sec. (final curing 16 min.) |
| - Spektramat (Ivoclar)                       | ca. 60 sec. (final curing 20 min.) |
| - Lampadaplust with light 71- 86W (Micerium) | ca. 10 min. (final curing 30 min.) |

Dental office curing times:

- |                                  |             |
|----------------------------------|-------------|
| - Translux CL (Kulzer)           | ca. 40 sec. |
| - Nou-Lite halogen lamp (Nouvag) | ca. 40 sec. |
| - ENA Lux Led (Micerium)         | ca. 20 sec. |

## USE AND STORAGE

Do not store above 25°C.

Do not use the product after the expiration date (see label on syringe or on "tips").

Due to hygienic reasons ENA HRi "Tips" and flow application needles should be used only once.

Use the material at room temperature. Medical device, for dental use only: keep away from children.

To avoid material waste, turn back the spindle after removing the material.

After use, close container with cap and keep it closed. Avoid direct exposure to sunlight.

If the material is not completely cured, it may discolor, mechanical properties deteriorate and pulpal inflammation can occur.

## Veneers



*Adult central with old restorations*



*Finished preparation*



*After veneer has been seated*

## Onlays and Crowns



*Patient with fractured centrals*




*Finished onlay*



*Completed restorations*

## Posterior Crowns and Onlays



*Restorations fabricated using 2 dentin shades and 1 enamel shade*



## Implants and Combination Cases

