## Instructions for Use

1-800-667-9622

Synca Marketing Inc.  
337 Marion, Le Gardeur,  
QC., Canada, J5Z 4W8

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed dentist or dental technicians.

### Description

FIBER FORCE is a series of glass fiber pre-impregnated with light-curing methacrylate-resin (prepreg). FIBER FORCE can be custom shaped and then polymerized by light curing. Designed for a dental office or a dental lab, easy-to-use FIBER FORCE is incorporated into a removable denture to provide a solid and esthetic reinforcement. FIBER FORCE is compatible with all types of methacrylate resins (self-cure, heat-cure, pressure-cure and microwavable).

### Indications

FIBER FORCE is used for the structural reinforcement in the production or repair of removable dentures made from methacrylate resin. Different configurations and associated resins are available.

- FIBER FORCE configurations are selected and used according to the application and space available as indicated in the table below.
- A corresponding white or pink resin is used in conjunction with the selected configuration(s).

| Selection of the FIBER FORCE | Production of removable dentures (injected/pressed) | Repair and reinforcement of removable dentures | Extension and reinforcement of partial dentures |
|------------------------------|---|--|---|
| Pink Mesh<br>118 x 53mm      | ***   |  |   |
| Pink Braid<br>Ø 1mm          |   | ***<br><b>X 3</b>                              | ***<br><b>X 2</b>                               |
| Pink Braid<br>Ø 2mm          |   | ***  | ***   |
| Pink UD Fibers<br>Ø 2mm      |   | ***  | ***   |

\*\*\* : Advised use ("X" indicates number of FIBER FORCE to use)

### Contraindications

Allergies to methacrylates

### Side effects

With the current state of knowledge, there are no known side effects.

### Precautions

- Always wear gloves to avoid handling FIBER FORCE with fingers: non-cured resin may irritate fingers and eyes. Sensitization from long-term exposure to the resin should not be ruled out as a possible hazard.
- Non-cured FIBER FORCE should not come into contact with mucous tissues.
- Fabrication of FIBER FORCE must be done extra-orally.
- Follow instructions in order to completely cure FIBER FORCE.
- When cutting or grinding cured splints, wear protective glasses, masks, gloves, and use appropriate air evacuation: fiberglass particles may irritate skin. In case of irritation, discomfort may be relieved with mild soap and water (light rubbing).
- Once package is opened, FIBER FORCE must be used quickly without exposure to light (sunlight as example): risk of premature curing.

### Composition (% by weight)

|                                   | preshaped fibers | resin in syringe |
|-----------------------------------|------------------|------------------|
| Glass fiber                       | 45-55%           | 0%               |
| Urethane dimethacrylate           | 30-40%           | 70%-85%          |
| Triethylene glycol dimethacrylate | 3-4%             | 8%-12%           |
| Inorganic elements and pigments   | 4-5%             | 8%-12%           |
| Catalysts/stabilizers             | <1%              | < 1%             |

### FIBER FORCE implementation

#### A. Production of a methyl methacrylate removable denture: injection technique

- 1) It is necessary to anchor the FIBER FORCE mesh to the plaster model so that the mesh will not move during the injection process. To accomplish this, make 2 small anchor holes in the center of the plaster model (2 mm deep).
- 2) Apply a wax separator to the plaster model and the anchor holes. Allow to dry completely.
- 3) Heat and adapt a wax sheet (0.2mm – 0.6 mm) to the model. Create 8 holes in the wax along the gingival crest and base of the model. These holes will be used to create cleats that are necessary for the proper spacing of the FIBER FORCE mesh during the injection process.
- 4) Preventing exposure to light, overfill the 8 holes in the wax with FIBER FORCE pink resin.
- 5) Preventing exposure to light, remove the FIBER FORCE mesh from its blister package. Cut to approximate size with the protective film in place. Remove the protective film.
- 6) Place the FIBER FORCE mesh over the wax on the model. Ensure that it is properly centered.
- 7) Place the model in the SPLINTVAC and adjust the height as indicated in the SPLINTVAC instructions. As per SPLINTVAC instructions, close the SPLINTVAC and use the vacuum pump to adapt the FIBER FORCE mesh to the model.
- 8) Place the SPLINTVAC and model into a light curing unit and light cure FIBER FORCE. Follow curing times listed below.
- 9) Remove the model from the SPLINTVAC. Remove the cured FIBER FORCE mesh from the model and remove any remaining wax from the model and the mesh (boiling water/steam). Let dry completely.
- 10) Ensure that any spaces in the mesh are free of resin. If necessary, pierce openings with a hand tool.
- 11) As necessary, cut any unadapted parts of the mesh. For best esthetics, the mesh should be cut just before the intersection of the teeth – gingival margin. Use short, sturdy scissors. Set the mesh aside.
- 12) Pick up the model and apply a coat of wax separator using your usual technique. Let it dry completely.
- 13) Apply wax to the model using your usual technique. Place 2 injection sites next to the anchor holes. Continue flasking procedure as usual.
- 14) Before closing the flask, overfill in the anchor holes made at A1) with Pink FIBER FORCE resin to ensure an excess.
- 15) Place the FIBER FORCE mesh over the plaster model and in contact with the resin in the anchor holes; this must locally coat the mesh. As necessary add a few drops of resin over the mesh.
- 16) Light cure the FIBER FORCE resin in a light curing unit, respecting the times indicated below. The mesh is now anchored to the plaster model.
- 17) Start the injection process with your usual methyl methacrylate resin, as per manufacturer's instructions.
- 18) Finish the denture as per your usual process. The FIBER FORCE mesh must remain completely within the resin.

#### B. Production of a methyl methacrylate removable denture: press technique

- 1) Apply a wax separator to the plaster model. Allow to dry completely.
- 2) Heat and adapt a wax sheet (0.2mm – 0.6 mm) to the model. Create 6 holes in the wax along the gingival crest and base of the model. These holes will be used to create cleats that are necessary for the proper spacing of the FIBER FORCE mesh during the injection process.
- 3) Preventing exposure to light, overfill the 6 holes in the wax with FIBER FORCE Pink resin.
- 4) Preventing exposure to light, remove the FIBER FORCE mesh from its blister package. Cut to approximate size with the protective film in place. Remove the protective film.
- 5) Place the FIBER FORCE mesh over the wax on the model. Ensure that it is properly centered.
- 6) Place the model in the SPLINTVAC and adjust the height as indicated in the SPLINTVAC instructions. As per SPLINTVAC instructions, close the SPLINTVAC and use the vacuum pump to adapt the FIBER FORCE mesh to the model.
- 7) Place the SPLINTVAC and model into a light curing unit and light cure FIBER FORCE. Follow curing times listed below.
- 8) Remove the model from the SPLINTVAC. Remove the cured FIBER FORCE mesh from the model and remove any remaining wax from the model and the mesh (boiling water/steam). Let dry completely.
- 9) Pierce openings with a hand tool to ensure that any spaces in the mesh are free of resin.
- 10) As necessary, cut any unadapted parts of the mesh. For best esthetics, the mesh should be cut just before the intersection of the teeth – gingival margin. Use short, sturdy scissors. Set the mesh aside.
- 11) Pick up the model and apply a coat of wax separator. Let it dry completely.
- 12) Use your usual waxing and flasking procedure.
- 13) Before closing the flask, place the FIBER FORCE mesh over the plaster model. Place your usual methyl methacrylate resin and press as per manufacturer's instructions.
- 14) Finish the denture as per your usual process. The FIBER FORCE mesh must remain completely within the resin.

#### C. Repair and reinforcement of removable dentures

- 1) Inspect the broken denture. If necessary, remove sharp edges.
- 2) Repair the denture following your usual procedure.
- 3) Create a groove along the inside arch of the denture at a 30° to 90° angle from the fracture, parallel to the teeth line (use a 1,5 to 2 mm fissure bur).
- 4) Place FIBER FORCE pink resin along the bottom of the groove.
- 5) Preventing exposure to light, remove a braid or UD fiber FIBER FORCE from its blister pack, and use scissors to cut it to size while it is still in the protective film. Remove the protective film.
- 6) Place the braid or UD fiber along the base of the groove, ensuring that it is well incorporated into the resin. Several fibers can be incorporated to maximize reinforcement.
- 7) Completely fill the groove with FIBER FORCE pink resin, ensuring that the fibers are completely covered. Overfill the groove.
- 8) Light cure the FIBER FORCE resin in a light curing unit or by a handheld curing light. Ensure that curing times are respected (see below).
- 9) Finish and polish the denture using your usual technique. The FIBER FORCE fibers must be completely covered by resin.

**Note:** Any self-cure methacrylate resin may be used instead of FIBER FORCE resin. Place the braid or UD FIBER FORCE fiber along the bottom of the groove. Polymerize it in a light-curing unit or using a hand-held curing light. Remove the hardened FIBER FORCE and set aside. Apply the self-cure methacrylate resin along the bottom of the groove. Place the hardened FIBER FORCE fiber into the groove. Cover the fiber with more self-cure methacrylate resin. Cure according to manufacturer's directions.

#### D. Reinforcing an extension of a partial denture.

- 1) Create a groove along the inside arch of the denture, parallel to the teeth line, at the area proposed for the extension (use a 1,5 to 2 mm fissure bur, 1 to 2mm deep).
- 2) Place FIBER FORCE pink resin along the bottom of the groove.
- 3) Preventing exposure to light, remove a braid or UD FIBER FORCE fiber from its blister package, and use scissors to cut it to size while it is still in the protective film. Remove the protective film.
- 4) Place the braid or UD fiber along the base of the groove, ensuring that it is well incorporated into the resin. Create a bow at the area where the extension will be made. This bow will solidify the area where the tooth will be added. Several fibers can be incorporated to maximize reinforcement.

- 5) Completely fill the groove with FIBER FORCE pink resin, ensuring that the fibers are completely covered. Overfill the groove.
- 6) Light cure the FIBER FORCE resin in a light curing unit or by a handheld curing light. Ensure that curing times are respected (see below).
- 7) Cover the bow with FIBER FORCE resin; in the same operation, place and affix the tooth to the resin.
- 8) Light cure the FIBER FORCE resin in a light curing unit or by a handheld curing light. Ensure that curing times are respected (see below).
- 9) Finish and polish the denture using your usual technique. The FIBER FORCE fibers must be completely covered by resin.

| Type of lampe | LED 5W light | Halogen, 1100mW/cm <sup>2</sup> | Halogen, 550mW/cm <sup>2</sup> | Xenon strobe light, 250mW/cm <sup>2</sup> | Neon, 6800mW/cm <sup>2</sup> | Mercury Vapor (Arc) |
|---------------|--------------|---------------------------------|--------------------------------|---|------------------------------|---------------------|
| Required time | 30 seconds   | 40 seconds                      | 2 minutes                      | 4 minutes                                 | 10 minutes                   | 20 minutes          |

### Note

Product reserved exclusively for dental usage.

Keep out of the reach of children.

Single use product.

Prevent exposure to light.

Store between 12°C and 30°C, in its original closed packaging.

Attention: see instructions for use.

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1-800-667-9622

Synca Marketing Inc.  
337 Marion, Le Gardeur,  
QC., Canada, J5Z 4W8

Caution: U.S. Federal law restricts this device to sale by or on the order of a licensed dentist or dental technicians.

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