



## **TempSpan™ Clear Matrix Material**

**Q. What are the working and setting times of *TempSpan* Clear Matrix Material?**

**A.** *TempSpan* Clear Matrix material offers a 1 minute working time and 2 minute oral set time.

**Q. Why did Pentron Clinical create *TempSpan* Clear Matrix Material?**

**A.** *TempSpan* Clear Matrix Material is a clear pre-operative material that reproduces fine detail for the fabrication of provisional restorations. Therefore, you can take your preliminary impression with this clear material and use the matrix to fabricate a quick and reliable restoration. This material is clear to provide feasibility to light cure dual-cure provisional materials intraorally or extraorally.

**Q. Can I use *TempSpan* Clear Matrix Material as a general impression material?**

**A.** Yes, *TempSpan* Clear Matrix can be used as a monophasic material for all general impressions. However, we do not recommend that you use *TempSpan* Clear Matrix as a bite registration material.

**Q. Which impression gun do I use with *TempSpan* Clear Matrix Material?**

**A.** *TempSpan* Clear Matrix was designed to fit the DS50 1:1/2:1 impression dispenser gun.

**Q. Does this material help to reduce the oxygen inhibited layer associated with dual-cure composite materials?**

**A.** Yes, light-curing dual-cure provisional materials within a clear matrix minimizes and may eliminate the oxygen inhibited smear layer.

**Q. Is *TempSpan* Clear Matrix Material hydrophilic?**

**A.** Yes.

**Q. What kind of impression tray should I use with your *TempSpan* Clear Matrix Material?**

**A.** Pentron Clinical recommends the use of clear plastic trays. However, you may also use a triple tray.