

INSTRUCTIONS FOR USE:

Bond-1 Technique Instructions for Direct Bonding to tooth structure:

1. Prepare and clean preparation site.
2. Etch enamel and dentin with 37% phosphoric acid etchant for 20 seconds. Rinse thoroughly. Remove excess moisture by blotting or with a light stream of air for 2 seconds. Leave surface moist.
3. Dispense *Bond-1* Primer/Adhesive and apply two consecutive coats to tooth within 10 seconds, to ensure a shiny surface.
4. **Critical Step:** Remove solvent from *Bond-1* Primer/Adhesive without disturbing resin surface. We recommend a gentle stream of air for a minimum of ten (10) seconds. Hold air syringe 1 inch from site, positioned so as not to disturb resin surface. If surface shine disappears after drying, re-coat with *Bond-1* Primer/Adhesive and re-dry as above.
5. Light cure for 10 seconds.
6. For direct composite, re-apply *Bond-1* remaining on the brush tip to tooth. Air dry for 10 seconds (no light cure is necessary). Apply composite.

Permanently Cementing Metal Restorations

Materials needed: etchant, *Bond-1* Primer/Adhesive, Cement-It® Universal C&B™ Cement

TOOTH

1. Prepare and clean restoration site per usual technique.
2. Total etch enamel and dentin for 20 seconds. Rinse thoroughly. Blot dry or use a light two second air stream. **Important:** Leave surface moist. If tooth is dry, poor bond strength and postoperative sensitivity may result.
3. Mix *Bond-1* Dual Cure Activator with *Bond-1* Primer/Adhesive. Using a fully saturated brush tip each time, apply two coats of *Bond-1* Primer/Adhesive to tooth within 10 seconds to ensure a resin rich, shiny surface.
4. **Critical Step:** Remove solvent from *Bond-1* Primer/Adhesive without disturbing resin surface. We recommend a gentle stream of air for a minimum of ten (10) seconds. Hold air syringe 1 inch from site, positioned so as not to disturb resin surface. (Avoid excess *Bond-1* Primer/Adhesive in internal line angles or point angles.) If surface shine disappears after drying, re-coat with *Bond-1* Primer/Adhesive and re-dry as above.
5. Light cure for ten seconds. Resultant layer (approx. 8µm) will not affect the fit of indirect restorations.

METAL

1. Sandblast inside of restoration and clean thoroughly.
2. Apply one coat of *Bond-1* Primer/Adhesive onto sandblasted internal surface of restoration.
3. Dry thoroughly with a gentle stream of air for about 10 seconds. The removal of moisture and *Bond-1* Primer/ Adhesive solvent is essential to achieving high bond strength. Do not light cure.
4. Mix *Cement-It Universal C&B* Cement (or other self/dual-cure resin cement) and place into restoration.
5. Place restoration as usual and follow cement kit instructions.

Bonding Metal Posts and Core Buildups

Materials needed: etchant, *Bond-1* Primer/Adhesive, *Cement-It Universal C&B* Cement, Build-It® FR™ Fiber Reinforced Core Build-Up Material

POST CANAL

1. Prepare and clean canal per usual technique.
2. Total etch all enamel and dentin (including post canal) for 20 seconds. Rinse thoroughly. Remove excess water by blotting with a lint-free, dry 2x2 gauze or use a light two second stream of air. Remove excess water in the canal space by blotting with a paper point. **Important:** Leave surface moist. If tooth is dry, poor bond strength and postoperative sensitivity may result.
3. a. Mix one drop of *Bond-1* Dual Cure Activator with two drops of *Bond-1* Primer/Adhesive. b.) Flood the post canal with *Bond-1* Primer/Adhesive (See technique tip). c.) Using a fully saturated brush tip, apply *Bond-1* Primer/ Adhesive to etched dentin and enamel surfaces. d.) Remove excess *Bond-1* in root canal with a paper point.
4. **Critical Step:** Remove solvent from *Bond-1* Primer/Adhesive without disturbing resin surface. We recommend a gentle stream of air for a minimum of ten (10) seconds. Hold air syringe 1 inch from site, positioned so as not to disturb resin surface. If surface shine disappears after drying, re-coat with *Bond-1* Primer/Adhesive and re-dry as above.
5. Light cure for ten (10) seconds for the surface enamel and/or dentin.

POST

1. Sandblast post and clean thoroughly.
2. Apply one coat of *Bond-1* Primer/Adhesive onto metal post (in case of cast post and core apply *Bond-1* Primer/ Adhesive to internal surface of core as well. Remove excess from line angles) Blow thin and dry to remove solvent.
3. Light cure for ten seconds. **Important:** For tight fitting posts and/or cores, light curing is not necessary.
4. Mix equal amounts of *Cement-It Universal C&B* Cement (or other self-cure resin cement) and place onto post and/or cast core.
5. Place post as usual and follow cement kit instructions.
6. If necessary, build or syringe composite core material.

For core buildups only

1. Thoroughly remove any temporary cement.
2. Follow instructions for direct composite, applying *Bond-1* Primer/Adhesive to the metal post as well.



INSTRUCTIONS FOR USE:

Permanently Cementing Porcelain and Indirect Composite Restorations

Materials needed: etchant, *Bond-1* Primer/Adhesive, Silane, *Lute-It*® Dual-Cure Cement

TOOTH

1. Remove temporary restoration. Prepare and clean restoration site per usual technique.
2. Try-in restoration for fit, color, etc. (See instructions in *Lute-It* Dual cure Cementation Kit.)
3. Total etch enamel and dentin for 20 seconds. Rinse thoroughly. Remove excess water by blotting or a two second stream of air. Important: Leave surface moist. If tooth is dry, poor bond strength and sensitivity may result.
4. Mix *Bond-1* Dual Cure Activator with *Bond-1* Primer/Adhesive. Using a fully saturated brush tip each time, apply two coats of *Bond-1* Primer/Adhesive to tooth within 10 seconds to ensure a resin rich, shiny surface.
5. Critical Step: Remove solvent from *Bond-1* Primer/Adhesive without disturbing resin surface. We recommend a gentle stream of air for a minimum of ten (10) seconds. Hold air syringe 1 inch from site, positioned so as not to disturb resin surface. (Avoid excess *Bond-1* Primer/Adhesive in internal line angles or point angles.)

PORCELAIN OR INDIRECT COMPOSITE

1. Apply Silane to inside of etched or sandblasted restoration. Allow Silane to stand for 60 seconds, then blow dry. (For composite restorations, the application of silane is optional.)
2. Brush on 1 coat of *Bond-1* Primer/Adhesive.
3. Remove solvent from *Bond-1* Primer/Adhesive with a jet of air. Light curing is not necessary.
4. Mix and apply proper shade of *Lute-It* Dual-cure Cement (other resin cement may be used) to inside of restoration.
5. Seat restoration and cure following manufacturer's technique instructions.

REPAIR

Composite to Metal

Materials needed: intraoral sandblaster, hydrofluoric acid gel, *Lute-It* Silane, *Bond-1* Primer/Adhesive, *Lute-It* Opaque Universal Stain.

1. If necessary, use a diamond bur to create a clean porcelain finish line. Rinse and dry.
2. Place rubber dam.
3. Sandblast site of broken porcelain and metal. Rinse and dry. If sandblasting is not possible, use a fine diamond bur to roughen up porcelain margin and exposed metal surface.
4. Brush hydrofluoric acid gel onto porcelain margin. Wait 1 minute. Rinse well and dry.
5. Apply *Lute-It* Silane to porcelain. Dry with jet of air.
6. Apply 1 coat of *Bond-1* Primer/Adhesive onto both the treated porcelain and the sandblasted metal. Dry thoroughly with a (minimum 10 seconds) jet of air.
7. Light cure *Bond-1* Primer/Adhesive for 10 seconds.
8. Apply *Lute-It* Opaque Universal Stain in a thin, even layer to metal surface. Light cure for 40 seconds.
9. Build and cure composite following manufacturer's instructions.

Composite to Composite

Materials needed: *Bond-1* Primer/Adhesive, intraoral sandblaster

1. Using fine diamond bur or carbide, grind existing composite to expose fresh surface.
2. Apply rubber dam and roughen the surface with an intraoral sandblaster, if possible. Rinse and dry with jet of air.
3. Apply 1 coat of *Bond-1* Primer/Adhesive to composite surface. Dry thoroughly.
4. Light cure for 10 seconds. Build and cure composite following manufacturer's instructions.

Composite to Porcelain

Materials needed: intraoral sandblaster, hydrofluoric acid etchant, silane, *Bond-1* Primer/Adhesive

1. Using a fine diamond bur or carbide, create a finish line in the porcelain.
2. Apply rubber dam and blast with intraoral sandblaster, if possible. Rinse and dry with jet of air.
3. Brush hydrofluoric acid gel onto porcelain margin. Wait 1 minute. Rinse well and dry.
4. Apply silane to porcelain. Air dry.
5. Apply 1 or 2 coats of *Bond-1* Primer/Adhesive to porcelain surface. Dry thoroughly.
6. Light cure for 10 seconds. Build and cure composite following manufacturer's instructions.

Using *Bond-1* Dual Cure Activator

Materials needed: etchant, *Bond-1* Primer/Adhesive, *Lute-It* or *Cement-It Universal C&B*, *Build-It F.R.*, direct composite

Bond-1 Dual Cure Activator is designed for use in conjunction with *Bond-1* and *Bond-1* C&B™ Primer/Adhesive. It is for use underneath a self-cured or dual-cured luting cement, core build-up composite or tooth filling composite where light curing is not indicated. *Bond-1* Dual Cure Activator is also recommended for situations where light curing to polymerize the bonding agent is limited (in a root canal, for example).

1. Follow the routine bonding procedure with *Bond-1* Primer/Adhesive and prepare the bonding surface.
2. Dispense one drop of *Bond-1* Dual Cure Activator and two drops of *Bond-1* Primer/Adhesive into a mixing well and briefly mix with a microbrush tip and apply directly to the bonding surface, the same way *Bond-1* Primer/ Adhesive is applied.
3. Dry the surface with a jet of air to remove the solvent in the bonding agent. The bonding surface should appear shiny and resin rich. If not, reapply the resin mix and re-dry again, if necessary.
4. Apply a self-cured or dual-cured luting cement (*Lute-It* or *Cement-It Universal C&B* cement), core build-up composite (*Build-It F.R.* Core Build-up Material) directly onto the self-curable bonding agent surface and allow to set.
5. Optionally, the material may be light cured for 10 seconds, with a standard halogen lamp, to achieve instant cure.

Pentron, *Bond-1*, *Cement-It*, *Build-It* and *Lute-It* are registered trademarks of Pentron® Clinical Technologies, LLC. *Universal C&B*, *F.R.* and *Bond-1* C&B are registered trademarks of Pentron Clinical Technologies, LLC.

