



Fini™ Polishing Disks

Q. How many sizes are available for the *Fini* polishing disks?

A. *Fini* polishing disks come in 9 sizes, with the choice of 4 grits and several different disk diameters.

Q. What advantages do *Fini* polishing disks offer over traditional polishing disks?

A. *Fini* polishing disks are very unique in that they do not possess a metal or plastic center ring. This eliminates the gouging and discoloration that can result from other polishing systems. *Fini* disks also have a larger usable surface area and are extremely flexible to provide the clinician easy access to interproximal and marginal areas.

Q. At what rotational speed should I run my handpiece to achieve best results?

A. For best results, use *Fini* polishing disks between 10,000 and 15,000 rpm. Do not exceed 25,000 rpm, or separation of the disk from the collar may occur.