








## Impression Material Comparison Chart

Viscosity	Description	Indication	Product	Work Time	Intraoral Time	Total Set Time	Color
<b>Thick'N Thin® Light Body</b>	A unique, thixotropic light-viscosity material that initially stacks on the tooth when syringed and then flows into the sulcus to capture critical subgingival detail.	W	<b>Correct Plus® Material</b>	1:00	3:30	4:30	
			<b>Correct Plus® Fast Set Material</b>	0:45	1:45	2:30	
<b>Light Body</b>	A traditional light-viscosity wash material characterized by adequate flowability.	W	<b>Correct Plus® Material</b>	1:00	3:30	4:30	
<b>Universal™ Body</b>	Truly an all-purpose impression material. Use it as a tray material and/or light body. A perfect substitute for alginate for night guards, bleach trays, or any other impression procedure.	W & C	<b>Correct Plus® Material</b>	1:00	3:30	4:30	
			<b>Correct Plus® Fast Set Material</b>	0:45	1:45	2:30	
			<b>Correct Plus® Mighty Mix™ Material</b>	1:00	3:30	4:30	
			<b>Correct Plus® Fast Set Mighty Mix™ Material</b>	0:45	1:45	2:30	
<b>Tray Material</b>	A heavy-viscosity material with a mousse-like consistency and high-end hardness. Ideal for Three-in-One impression trays.	Q, FA, TIO	<b>Correct Plus® Material (Unflavored/Berry Flavored)</b>	1:00	3:30	4:30	
			<b>Correct Plus® Mighty Mix™ Material</b>	1:00	3:30	4:30	
<b>Auto-Mix Putty</b>	A heavy viscosity material that provides the hydraulic drive of a hand-mix putty in a convenient auto-mix cartridge.	Q, FA, TIO	<b>Correct Plus® Material (Unflavored/Berry Flavored)</b>	1:00	3:30	4:30	
			<b>Correct Plus® Fast Set Material</b>	0:45	1:45	2:30	
			<b>Correct Plus® Mighty Mix™ Material</b>	1:00	3:30	4:30	
			<b>Correct Plus® Fast Set Mighty Mix™ Material</b>	0:45	1:45	2:30	
<b>Putty</b>	An easy-to-mix, heavy-viscosity putty with a high-end hardness to resist distortion and provide maximum hydraulic drive of the light body material.	Q & FA	<b>Correct Plus® Material (Unflavored/Berry Flavored)</b>	1:00	3:30	4:30	
			<b>Correct VPS® Material</b>	1:00	3:30	4:30	

## Bite Registration Material Comparison Chart

Viscosity	Description	Indication	Product	Work Time	Intraoral Time	Total Set Time	Color
<b>Bite Material/ Bite Superfast</b>	These bite registration materials have a mousse-like viscosity and high-end hardness. When set, its rigidity facilitates trimming.	Bite Registration	<b>Correct Plus® Material</b>	0:15	1:00	1:15	
			<b>Correct Quick® Material</b>	0:30	1:00	1:30	

Indication Abbreviations: W=Wash, C=Custom Tray, Q=Quadrant Tray, FA=Full Arch Tray, TIO=Three-in-One Impression Tray